remel

BactiDrop™ Dobell & O'Connor Iodine

(English)

INTENDED USE

Remel BactiDrop™ Dobell & O'Connor Iodine stain is recommended for use to prepare wet mounts of intestinal protozoa.

SUMMARY AND EXPLANATION

This iodine solution was formulated by Dobell and O'Connor in 1921 for the staining of protozoan cysts in fecal specimens.¹ They found that a weak iodine solution did not coagulate fecal particles that could obscure or entrap organisms. It has also been reported that a weak iodine solution does not adversely affect the refractile nature of stained organisms, such as concentrated iodine does.²

PRINCIPLE

The iodine mount is useful for the examination of cysts and eggs, but trophozoites are killed. The chromatin material of amebic cysts stain dark brown against the yellow-brown cytoplasm, the nuclear structures are differentiated and the glycogen masses stain a mahogany brown.²

REAGENTS (CLASSICAL FORMULA)*

Potassium Iodide (CAS 7681-11-0) 20.0 g
Iodine (CAS 7553-56-2) 10.0 g
Demineralized Water (CAS 7732-18-5) 1000.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

CAUTION! May cause allergic skin reaction. May cause eye, skin, and respiratory tract irritation.

This product is For *In Vitro* Diagnostic Use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after use. Directions should be read and followed carefully. Refer to Material Safety Data Sheet for additional information.

STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 20-25°C until used. Do not freeze or overheat. Protect from light.

PRODUCT DETERIORATION

This product should not be used if (1) there is evidence of dehydration, (2) the color has changed, (3) the expiration date has passed, or (4) there are other signs of deterioration. The expiration date applies to the product in its intact container when stored as directed. Discard remaining portion of partially used ampule at end of workday.

SPECIMEN COLLECTION, STORAGE, TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{3,4}

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Applicator sticks, plain and cotton-tipped, (4) Incubator, slide warmer, (5) Quality control organisms, (6) Glass slides, coverslips, (7) Microscope with a calibrated ocular micrometer, (8) Disposable glass or plastic pipettes.

PROCEDURE

Place dropper in the assembled, reusable ampule crusher provided. Hold the dropper/crusher in an upright position and lightly tap the bottom to dislodge any bubbles that may have formed. Grasp the middle of the dropper/crusher with the thumb and forefinger, and with the tip pointing away, press gently to crush the ampule. Invert dropper and squeeze slightly to dispense in a dropwise fashion.

Test Procedure:

Dispense stain as needed for procedures involving examination of intestinal protozoa.

QUALITY CONTROL

All lot numbers of BactiDrop™ Dobell & O'Connor Iodine have been tested using the following quality control organisms and have been found to be acceptable. Testing with a known, fresh, reference specimen should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

CONTROL

Giardia intestinalis
ATCC® 30888

RESULTS

Intact protozoa, mahogany brown in color

LIMITATIONS

1. Iodine preparations fade quickly if exposed to light.³
2. The karyosome of *Iodamoeba butschlii* does not stain with iodine.
3. The glycogen of older cysts may not stain with iodine.²

4. Specimens collected from patients receiving barium sulfate, mineral oil, bismuth, nonabsorbable antidiarrheal preparations, antimalarials, and certain antibiotics are unsuitable for examination.²
5. A minimum of three fecal specimens is recommended for routine examination for parasites prior to treatment.²
6. Protozoan organisms should not be identified on the basis of wet mount (with few exceptions). Permanent stained smears should be examined to confirm the identification.²

BIBLIOGRAPHY

1. Dobell, C. and F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute (CLSI)/NCCLS. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract; Approved Guideline. 2nd ed. M28-A2. CLSI, Wayne, PA.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Tenover, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.
4. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11th ed. Mosby, Inc., St. Louis, MO.

PACKAGING

BactiDrop™ Dobell & O'Connor Iodine (0.75 ml/Ampule):
REF 21510..... 50 Ampules/Pk

Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

BactiDrop™ is a trademark of Remel Inc.
ATCC® is a registered trademark of American Type Culture Collection.
CAS (Chemical Abstracts Service Registry No.)

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BIBLIOGRAPHIE

1. Dobell, C. and F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute (CLSI)/NCCLS. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract. Approved Guideline, 2nd ed. M28-A2. CLSI, Wayne, PA.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.

4. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11th ed. Mosby, Inc., St. Louis, MO.

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BIBLIOGRAPHIE

1. Dobell, C. und F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute (CLSI)/NCCLS. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract. Approved Guideline, 2nd ed. M28-A2. CLSI, Wayne, PA.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Tenover und R.H. Tenover. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.

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BIBLIOGRAFIA

1. Dobell, C. and F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute (CLSI)/NCCLS. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract; Approved Guideline. 2nd ed. M28-A2. CLSI, Wayne, PA.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.

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1. Dobell, C. and F.W. O'Connor. 1921. Intestinal Protozoa of Man. William Wood, New York.
2. Clinical and Laboratory Standards Institute (CLSI)/NCCLS. 2005. Procedures for the Recovery and Identification of Parasites from the Intestinal Tract. Approved Guideline, 2nd ed. M28 A2. CLSI, Wayne, PA.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8th ed. ASM, Washington, D.C.

4. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11th ed. Mosby, Inc., St. Louis, MO.

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BactiDrop™ es una marca comercial de Remel Inc.